

Phospho-VCL (Y821) Antibody



PACO05292

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

VCL

Uniprot

P18206

Synonyms:

VCL; Vinculin; Metavinculin

Immunogen:

Synthesized peptide derived from human Vinculin around the phosphorylation site of Y821.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A