

VHL Antibody



PACO05298

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

VHL

Uniprot

P40337

Synonyms:

VHL; Von Hippel-Lindau disease tumor suppressor; Protein G7; pVHL

Immunogen:

Synthesized peptide derived from human VHL around the non-phosphorylation site of S68.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A