

Phospho-EZR (Y353) Antibody



PACO05299

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

EZR

Uniprot

P15311

Synonyms:

EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81

Immunogen:

Synthesized peptide derived from human Ezrin around the phosphorylation site of Y353.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A