

Phospho-VIM (Y61) Antibody



PACO05301

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

VIM

Uniprot

P08670

Synonyms:

VIM; Vimentin

Immunogen:

Synthesized peptide derived from human Vimentin around the phosphorylation site of Y61.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A