

VIP Antibody



PACO05302

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:

Gene ID:

VIP

Uniprot

P01282

Synonyms:

VIP; VIP peptides

Immunogen:

Synthesized peptide derived from the Internal region of human VIP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A