

# ZAP70 Antibody



PACO05325

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ZAP70

**Uniprot**

P43403

**Synonyms:**

ZAP70; SRK; Tyrosine-protein kinase ZAP-70; 70 kDa zeta-chain associated protein; Syk-related tyrosine kinase

**Immunogen:**

Synthesized peptide derived from human ZAP-70 around the non-phosphorylation site of Y315.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A