

ZIC1/ZIC2/ZIC3/ZIC4/ZIC5 Antibody



PACO05329

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:40000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ZIC1/ZIC2/ZIC3/ZIC4/ZIC5

Uniprot

Q15915/O95409/O60481/Q8N9L1/Q96T25

Synonyms:

ZIC1; ZIC; ZNF201; Zinc finger protein ZIC 1; Zinc finger protein 201; Zinc finger protein of the cerebellum 1; ZIC2; Zinc finger protein ZIC 2; Zinc finger protein of the cerebellum 2; ZIC3; ZNF203; Zinc finger protein ZIC 3; Zinc finger p

Immunogen:

Synthesized peptide derived from the Internal region of human Zic1/2/3/4/5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A