

CA1 Antibody



PACO05332

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

CA1

Uniprot

P00915

Synonyms:

CA1; Carbonic anhydrase 1; Carbonate dehydratase I; Carbonic anhydrase B; CAB;
Carbonic anhydrase I; CA-I

Immunogen:

Synthesized peptide derived from the Internal region of human CA I.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A