

# Phospho-ZYX (S142) Antibody



PACO05342

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

ZYX

**Uniprot**

Q15942

**Synonyms:**

ZYX; Zyxin; Zyxin-2

**Immunogen:**

Synthesized peptide derived from human Zyxin around the phosphorylation site of S142.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A