

**Phospho-
TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA3E/TUB
A4A (Y272) Antibody
PACO05344**



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA3E/TUBA4A

Uniprot

Q71U36/P68363/Q9BQE3/Q13748/Q6PEY2/P68366

Synonyms:

TUBA1A; TUBA3; Tubulin alpha-1A chain; alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; alpha-tubulin ubiquitous; Tubulin K-alpha-1; Tubulin alpha-ubiquitous chain; TUBA1C; TUBA6; Tubulin alpha-1C c

Immunogen:

Synthesized peptide derived from human Tubulin α 1/3/4 around the phosphorylation site of Y272.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A