

FZD5/FZD8 Antibody



PACO05347

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

FZD5/FZD8

Uniprot

Q13467/Q9H461

Synonyms:

FZD5; C2orf31; Frizzled-5; Fz-5; hFz5; FzE5; FZD8; Frizzled-8; Fz-8; hFz8

Immunogen:

Synthesized peptide derived from the N-terminal region of human Frizzled-5/8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A