

Phospho-CAD (T456) Antibody



PACO05349

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

CAD

Uniprot

P27708

Synonyms:

CAD; CAD protein

Immunogen:

Synthesized peptide derived from human CAD around the phosphorylation site of T456.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A