

FKRP Antibody



PACO05353

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

FKRP

Uniprot

Q9H9S5

Synonyms:

FKRP; Fukutin-related protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human FKRP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A