

OR5H6 Antibody



PACO05357

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5H6

Uniprot

Q8NGV6

Synonyms:

OR5H6; Olfactory receptor 5H6; Olfactory receptor OR3-11

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5H6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A