

# OR51G1 Antibody



PACO05358

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

OR51G1

**Uniprot**

Q8NGK1

**Synonyms:**

OR51G1; OR51G3P; Olfactory receptor 51G1; Olfactory receptor 51G3; Olfactory receptor OR11-29

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 51G1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A