

CALB2 Antibody



PACO05360

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

CALB2

Uniprot

P22676

Synonyms:

CALB2; CAB29; Calretinin; CR; 29 kDa calbindin

Immunogen:

Synthesized peptide derived from the Internal region of human Calretinin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A