

EPS8L3 Antibody



PACO05363

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

EPS8L3

Uniprot

Q8TE67

Synonyms:

EPS8L3; EPS8R3; Epidermal growth factor receptor kinase substrate 8-like protein 3; EPS8-like protein 3; Epidermal growth factor receptor pathway substrate 8-related protein 3; EPS8-related protein 3

Immunogen:

Synthesized peptide derived from the Internal region of human Eps8L3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A