

ABHD8 Antibody



PACO05364

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

ABHD8

Uniprot

Q96I13

Synonyms:

ABHD8; Abhydrolase domain-containing protein 8

Immunogen:

Synthesized peptide derived from the Internal region of human ABHD8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A