

CALCA/CALCA Antibody



PACO05367

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CALCA/CALCA

Uniprot

P06881/P01258

Synonyms:

CALCA; CALC1; Calcitonin gene-related peptide 1; alpha-type CGRP; Calcitonin gene-related peptide I; CGRP-I; CALCA; CALC1; Calcitonin

Immunogen:

Synthesized peptide derived from the Internal region of human Calcitonin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A