

# ZNF668 Antibody



PACO05374

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

ZNF668

**Uniprot**

Q96K58

**Synonyms:**

ZNF668; Zinc finger protein 668

**Immunogen:**

Synthesized peptide derived from the Internal region of human ZNF668.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A