

MTERF1 Antibody



PACO05376

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

MTERF1

Uniprot

Q99551

Synonyms:

MTERF; Transcription termination factor; mitochondrial; Mitochondrial transcription termination factor 1; mTERF

Immunogen:

Synthesized peptide derived from the Internal region of human mTERF.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A