

MAP9 Antibody



PACO05381

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

MAP9

Uniprot

Q49MG5

Synonyms:

MAP9; ASAP; Microtubule-associated protein 9; Aster-associated protein

Immunogen:

Synthesized peptide derived from the Internal region of human MAP-9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A