

CEP76 Antibody



PACO05385

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

CEP76

Uniprot

Q8TAP6

Synonyms:

CEP76; C18orf9; Centrosomal protein of 76 kDa; Cep76

Immunogen:

Synthesized peptide derived from the Internal region of human CEP76.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A