

PACO05397

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

CUBN

**Uniprot**

O60494

**Synonyms:**

CUBN; IFCR; Cubilin; 460 kDa receptor; Intestinal intrinsic factor receptor; Intrinsic factor-cobalamin receptor; Intrinsic factor-vitamin B12 receptor

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Cubilin.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A