

APOL5 Antibody



PACO05406

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

APOL5

Uniprot

Q9BWW9

Synonyms:

APOL5; Apolipoprotein L5; Apolipoprotein L-V; ApoL-V

Immunogen:

Synthesized peptide derived from the C-terminal region of human ApoL5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A