

OR11H1/OR11H2/OR11H12 Antibody



PACO05410

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR11H1/OR11H2/OR11H12

Uniprot

Q8NG94/Q8NH07/B2RN74/Q6IEX0

Synonyms:

OR11H1; Olfactory receptor 11H1; Olfactory receptor OR22-1; OR11H2; C14orf15; OR11H2P; Olfactory receptor 11H2; Olfactory receptor OR14-1; OR11H12; Olfactory receptor 11H12Q6IEX0

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 11H1/11H2/11H12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A