

OR51E2 Antibody



PACO05414

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR51E2

Uniprot

Q9H255

Synonyms:

OR51E2; PSGR; Olfactory receptor 51E2; HPRAJ; Olfactory receptor OR11-16; Prostate-specific G-protein coupled receptor

Immunogen:

Synthesized peptide derived from the C-terminal region of human PSGR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A