

# OR2C3 Antibody



PACO05420

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

OR2C3

**Uniprot**

Q8N628

**Synonyms:**

OR2C3; OR2C4; OR2C5P; Olfactory receptor 2C3; Olfactory receptor 2C4; Olfactory receptor 2C5; Olfactory receptor OR1-30

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2C3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A