

GPR63 Antibody



PACO05421

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:20000, IF:1:200-1:1000

Protein Background:**Gene ID:**

GPR63

Uniprot

Q9BZJ6

Synonyms:

GPR63; PSP24B; Probable G-protein coupled receptor 63; PSP24-2; PSP24-beta

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR63.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A