

CAMK2A/CAMK2B/CAMK2D Antibody



PACO05422

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CAMK2A/CAMK2B/CAMK2D

Uniprot

Q9UQM7/Q13554/Q13557

Synonyms:

CAMK2A; CAMKA; KIAA0968; Calcium/calmodulin-dependent protein kinase type II subunit alpha; CaM kinase II subunit alpha; CaMK-II subunit alpha; CAMK2B; CAM2; CAMK2; CAMKB; Calcium/calmodulin-dependent protein kinase type II subunit beta; Ca

Immunogen:

Synthesized peptide derived from human CaMKII α ; / β ; / δ around the non-phosphorylation site of T305.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A