

STMN4 Antibody



PACO05423

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

STMN4

Uniprot

Q9H169

Synonyms:

STMN4; Stathmin-4; Stathmin-like protein B3; RB3

Immunogen:

Synthesized peptide derived from the C-terminal region of human RB3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A