

ADAMTS12 Antibody



PACO05429

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

ADAMTS12

Uniprot

P58397

Synonyms:

ADAMTS12; A disintegrin and metalloproteinase with thrombospondin motifs 12; ADAM-TS 12; ADAM-TS12; ADAMTS-12

Immunogen:

Synthesized peptide derived from the Internal region of human ADAMTS-12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A