

OR12D3 Antibody



PACO05430

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

OR12D3

Uniprot

Q9UGF7

Synonyms:

OR12D3; Olfactory receptor 12D3; Hs6M1-27; Olfactory receptor OR6-27

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 12D3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A