

MIA Antibody



PACO05432

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

MIA

Uniprot

Q16674

Synonyms:

MIA; Melanoma-derived growth regulatory protein; Melanoma inhibitory activity protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human MIA.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A