

FZD4 Antibody



PACO05441

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:5000, IF:1:200-1:1000

Protein Background:**Gene ID:**

FZD4

Uniprot

Q9ULV1

Synonyms:

FZD4; Frizzled-4; Fz-4; hFz4; FzE4; CD antigen CD344

Immunogen:

Synthesized peptide derived from the Internal region of human Frizzled-4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A