

PCDH11Y/PCDH11X Antibody



PACO05447

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

PCDH11Y/PCDH11X

Uniprot

Q9BZA8/Q9BZA7

Synonyms:

PCDH11Y; PCDH11; PCDH22; PCDHY; Protocadherin-11 Y-linked; Protocadherin-11; Protocadherin on the Y chromosome; PCDH-Y; Protocadherin prostate cancer; Protocadherin-PC; Protocadherin-22; PCDH11X; KIAA1326; PCDH11; PCDHX; Protocadherin-11 X-

Immunogen:

Synthesized peptide derived from the Internal region of human Protocadherin-11.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A