

Cleaved-CASP3 (D175) Antibody



PACO05465

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

CASP3

Uniprot

P42574

Synonyms:

CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1

Immunogen:

Synthesized peptide derived from the Internal region of human Caspase-3 p17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A