

RBP5 Antibody



PACO05476

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

RBP5

Uniprot

P82980

Synonyms:

RBP5; Retinol-binding protein 5; Cellular retinol-binding protein III; CRBP-III; HRBPiso

Immunogen:

Synthesized peptide derived from the N-terminal region of human CRBP-III.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A