

ZNF541 Antibody



PACO05506

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ZNF541

Uniprot

Q9H0D2

Synonyms:

ZNF541; Zinc finger protein 541

Immunogen:

Synthesized peptide derived from the C-terminal region of human ZNF541.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A