

PACO05513

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

RECK

Uniprot

O95980

Synonyms:

RECK; ST15; Reversion-inducing cysteine-rich protein with Kazal motifs; hRECK; Suppressor of tumorigenicity 15 protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human RECK.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A