

C2orf40 Antibody



PACO05519

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

C2orf40

Uniprot

Q9H1Z8

Synonyms:

C2orf40; ECRG4; Augurin; Esophageal cancer-related gene 4 protein

Immunogen:

Synthesized peptide derived from the Internal region of human ECRG4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A