

GPR123 Antibody



PACO05521

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:5000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ADGRA1

Uniprot

Q86SQ6

Synonyms:

GPR123; KIAA1828; Probable G-protein coupled receptor 123

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR123.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A