

TMEM185A Antibody



PACO05523

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

TMEM185A

Uniprot

Q8NFB2

Synonyms:

TMEM185A; CXorf13; FAM11A; Transmembrane protein 185A; Protein FAM11A

Immunogen:

Synthesized peptide derived from the C-terminal region of human TMEM185A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A