

# PAR6B Antibody



PACO05527

---

## Product Information

<b>Size:</b> 50ug	<b>Protein Background:</b>
<b>Reactivity:</b> Human, Mouse	<b>Gene ID:</b> PAR6B
<b>Source:</b> Rabbit	<b>Uniprot</b> Q9BYG5
<b>Isotype:</b> IgG	<b>Synonyms:</b> PAR6B; PAR6B; Partitioning defective 6 homolog beta; PAR-6 beta; PAR-6B
<b>Applications:</b> ELISA, WB, IHC	<b>Immunogen:</b> Synthesized peptide derived from the Internal region of human PAR-6 &beta; .
<b>Recommended dilutions:</b> ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300	<b>Storage:</b> Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A