

EMR3 Antibody



PACO05532

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

ADGRE3

Uniprot

Q9BY15

Synonyms:

EMR3; EGF-like module-containing mucin-like hormone receptor-like 3; EGF-like module receptor 3

Immunogen:

Synthesized peptide derived from the C-terminal region of human EMR3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A