

SFT2D3 Antibody



PACO05539

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

SFT2D3

Uniprot

Q58719

Synonyms:

SFT2D3; Vesicle transport protein SFT2C; SFT2 domain-containing protein 3

Immunogen:

Synthesized peptide derived from the Internal region of human SFT2D3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A