

# ABHD14B Antibody



PACO05540

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ABHD14B

**Uniprot**

Q96IU4

**Synonyms:**

ABHD14B; CIB; Abhydrolase domain-containing protein 14B; CCG1-interacting factor B

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human ABHD14B.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A