

# GALR3 Antibody



PACO05541

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

GALR3

**Uniprot**

O60755

**Synonyms:**

GALR3; GALNR3; Galanin receptor type 3; GAL3-R; GALR-3

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human GALR3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A