

RGS5 Antibody



PACO05542

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

RGS5

Uniprot

O15539

Synonyms:

RGS5; Regulator of G-protein signaling 5; RGS5

Immunogen:

Synthesized peptide derived from the Internal region of human RGS5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A