

LRP11 Antibody



PACO05544

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

LRP11

Uniprot

Q86VZ4

Synonyms:

LRP11; Low-density lipoprotein receptor-related protein 11; LRP-11

Immunogen:

Synthesized peptide derived from the Internal region of human LRP11.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A